



EZ Read 800 ELISA and Galapagos software for EZ Read delivers everything you need to make your absorbance assays (e.g. ELISA, Bradford and cell viability assays) easy.

VERSATILE DESIGN FOR RESEARCH LABORATORIES

Protein Assay	Filter	Included
Bradford	595nm	Yes
BCA	562nm	Yes
Cell Viability Assays		
MTT	492nm	Yes
Cell Density	620nm	Yes
ELISA Assay Substrate		
PNPP, ABTS	405nm	Yes
OPD	492nm	Yes
TMB	450nm	Yes
Reference wavelength	620nm	Yes

KEY FEATURES

Galapagos for EZ Read Software makes 'Click and drag' approach makes usability a key feature.

- Easy to create new methods.
- Select only the data you want for data reporting
- So easy you won't need to export your data into Excel.
- Windows 7 & 8 compatibility as well as older Windows OS from XP.

EZ Read 800 ELISA

- <5sec measurement of 96-well plate
- Measures **flat, round and v-bottomed 96-well** microplates with standard ANSI SBS footprints are measured
- **Accurate measurements** ensured by automatic lamp calibration prior to each reading
- Shaking for well homogeneity
- **USB 2.0 PC connection**

ORDERING INFORMATION

80-4002-01

EZ Read 800 Research Microplate Reader

SG010160

Quality Control Plate with WinQC software

Distributors Worldwide



77, avenue de Champagne
91420 MORANGIS

Tél. : 33 (0)1 60 19 69 29 Fax : 33 (0)1 60 19 69 25
www.bioserv.fr info@bioserv.fr

