# **EZ READ 800 Plus ELISA**





EZ Read 800 Plus ELISA and Galapagos software for EZ Read delivers everything you need to make your ELISAs easy.

### **DESIGNED FOR ELISA ASSAYS**

ELISA Assay Substrate	Filter	Included
PNPP	405nm	Yes
ABTS	405nm	Yes
OPD	492nm	Yes
Slow TMB	450nm	Yes
Turbo TMB	450nm	Yes
Ultra TMB	450nm	Yes
Reference wavelength	620nm	Yes

## **ONBOARD SOFTWARE**

- Thresholds, replicate verification and test validation formulae
- Stores 120 methods and the results of up to 100 plates
- Flexible control with user interface with multiple languages: French, German, Spanish and Portuguese

#### **KEY FEATURES**

# Galapagos for EZ Read Software's 'click and drag' approach makes usability a key feature.

- Easy to create new plate layouts
- Add multiple assays to a plate
- Add multiple plates per analysis
- Select only the data you want for data reporting
- Windows 7 & 8 compatibility as well as older Windows OS from XP.

#### EZ Read 800 Plus ELISA

- **<5sec** measurement of 96-well plate
- Measures flat, round and v-bottomed 96well microplates and stripwells with standard ANSI SBS footprints are measured
- Accurate measurements ensured by automatic lamp calibration prior to each reading
- Four **shaking** modes for complete mixing of solutions and more reproducible results.
- USB 2.0 PC connection

#### **ORDERING INFORMATION**

80-4002-03	EZ Read 800 ELISA Microplate Reader*
SG010160	Quality Control Plate with WinQC software

Distributors Worldwide



77, avenue de Champagne 91420 MORANGIS Tél. : 33 (0)1 60 19 69 29 Fax : 33 (0)1 60 19 69 25 www.bioserv.fr info@bioserv.fr

